

Code:
CCL2

Description:

20 Chemical Absorbent
Cushions, 30cm x 35cm -
Boxed

Absorbency / Capacity:
80 Litres

Dimensions:
30cm x 35cm

Application:
Chemical

Quantity
1

Weight - Kg:
9kg

Ean:
0617285799422



20 Chemical Absorbent Cushions, 30cm x 35cm - Boxed

20 Chemical absorbent cushions for use on acids, alkalis and strong chemicals. Ideal for use in drip trays to absorb leaked fluids easily and efficiently and to reduce further spillage during disposal of drip tray contents. Place directly under leaking hoses and fittings to absorb and contain liquids. Place on spilled liquids to absorb cleanly and efficiently.

Chemical

- Acids
- Alkalis
- Strong chemicals
- Solvents
- Not for use on oxidisers



T: 01487 823184
E: orders@fentex.co.uk
www.fentex.co.uk

